

## ***SO-LOW ECONOMY REFRIGERATOR***

### ***MODEL EC4-10GD***

**10 Cubic Ft. Upright – Temperature Range: 2°C to 8°C**

**Glass Door**

#### ***Specifications:***

* Capacity	10 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	24" x 26" x 54"
	60 x 66 x 137 (cm)
* Interior w x d x h	20.75" x 21" x 37.75"
	52 x 53 x 96 (cm)

#### ***Features:***

Temperature Control	Digital display
Access Port	1/2" Diameter Top Back Wall
Door	Single glass, locking
Interior finish	White coated steel
Legs	Four leveling
Shelves	Three Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Compressor	1/5 hp
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	2.0 / 15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	300 lbs.
Warranty-(Domestic)	Three year parts and labor
	Compressors additional
	2 years, no labor.

