

## ***SO-LOW ECONOMY REFRIGERATOR***

### ***MODEL EC4-12GD***

**12 Cubic Ft. Upright – Temperature Range: 2°C to 10°C**

**Glass Door**

#### ***Specifications:***

* Capacity	12 cubic ft.
* Temperature Range	2°C to 10°C
* Dimensions	
* Exterior w x d x h	24" x 24" x 62" 60 x 60 x 157 (cm)
* Interior w x d x h	20.75" x 21" x 45.75" 52 x 53 x 116 (cm)

#### ***Features:***

Door	Single glass, locking
Interior finish	White coated steel
Legs	Four leveling
Shelves	Four Adjustable
Interior Lighting	One
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/5 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0/15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	400 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice