

## ***SO-LOW ECONOMY REFRIGERATOR***

### ***MODEL EC4-27SD***

**27 Cubic Ft. (764 liters) Upright – Temperature Range: 2°C to 10°C**

**Solid Door**

#### ***Specifications:***

* Capacity	27 cubic ft. / 764 liters
* Temperature Range	2°C to 10°C
* Dimensions	
* Exterior w x d x h	30" x 34" x 78" 76 x 86 x 198 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

#### ***Features:***

Door	Single solid, locking
Interior finish	White Coated Steel
Legs	Four, leveling
Shelves	Five
Interior Lighting	One
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/3 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.8/15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice