

SO-LOW ECONOMY REFRIGERATOR MODEL EC4-38SGD

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

38 Cubic Ft. Upright – Temperature Range: 2°C to 8°C

Two Glass Sliding Doors

Specifications:

* Capacity	38 cubic ft. / 1076 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	43.5" x 32" x 78" 110 x 81 x 198 (cm)
* Interior w x d x h	40.5 x 27" x 55" 102 x 68 x 140 (cm)

Features:

Interior	White powder-coated Steel
Doors	Double glass – sliding Positive seal, self closing Lockable
Legs	Four, adjustable
Shelves	Ten
Interior Lighting	One
Insulation	Polyurethane
Refrigerant:	R134A
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps	8.8/15 Amp circuit required
Supply Plug	5-15 NEMA
Shipping Weight	550 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice