

SO-LOW ECONOMY REFRIGERATOR MODEL EC4-74GD

**72 Cubic Ft. (2095 liters) Upright – Temperature Range: 2°C to 8°C
Three Glass Hinged Doors**

Specifications:

* Capacity	72 cubic ft. / 2095 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	75" x 32" x 78" 190 x 81 x 198 (cm)
* Interior w x d x h	72" x 27" x 55" 183 x 69 x 140 (cm)

Features:

Doors	Triple glass, locking Positive seal, self closing
Shelves	Fifteen, adjustable
Legs	Four, leveling
Interior Lighting	Two
Insulation	Polyurethane CFC Free
Refrigerant	R134A CFC free
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	9.9/15
Supply Plug	NEMA 5-15P UL listed
Shipping Weight	850 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice