

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW UNDERCOUNTER REFRIGERATOR/ FREEZER COMBINATION MODEL MV18-6UCRF

6.1 Cubic Ft.(173 liters) Upright

Temperature Range: 1°C to 7°C (R)

-18°C to -9°C (F)

Specifications:

*Capacity	5.3 cubic ft. (R) .7 cubic ft. (F)
*Temperature Range	1°C to 7°C (R) -18°C to -9° (F)
*Dimensions	
*Exterior w x d x h	23-7/8" x 24" x 34" 61 x 61 x 86 (cm)

Features:

*COMMERCIAL UL LISTED

*MANUAL DEFROST CYCLE

*FORCED AIR VENTILATION
THROUGH FRONT GRILLE

*NO ADDITIONAL VENTILATION
CLEARANCE REQUIRED TO BUILD
IN UNDERCOUNTER

*R134A REFRIGERATION SYSTEM

*3 IN-DOOR SHELVES

*115 VOLT, 60 HERTZ, 1 PHASE

*3.3 AMPS

*WARRANTY, 1 YEAR P&L, ADDITIONAL
4 YEARS ON COMPRESSORS PARTS ONLY

*DOUBLE INSULATED POWER CORD WITH
THREE PRONG GREEN DOT HOSPITAL PLUG.

*EXTERIOR 22 GAUGE STEEL POWDER COAT
PAINT FINISH

*INTERIOR VACUUM FORM HIGH IMPACT POLYSTYRENE
SEAMLESS WITH ROUNDED CORNERS FOR EASY CLEANING



Specifications subject to change without notice