

## *SO-LOW UNDERCOUNTER FREEZER*

### *MODEL MV23-4UCF*

**4.5 Cubic Ft. (127 liters) Upright**

**Temperature Range: -23°C to -17°C**

**Solid Door**

#### *Specifications:*

* Capacity	4.5 cubic ft. / 127 liters
* Temperature Range	-23°C to -17°C
* Dimensions	
* Exterior w x d x h	23-7/8" x 24" x 34" 61 x 61 x 86 (cm)

#### *Features:*

- \*COMMERCIAL UL LISTED
- \*MANUAL DEFROST CYCLE
- \*FORCED AIR VENTILATION THROUGH FRONT GRILLE
- \*NO ADDITIONAL VENTILATION CLEARANCE REQUIRED TO BUILD IN UNDERCOUNTER
- \*R134A REFRIGERATION SYSTEM
- \*FLAT DOOR LINER
- \*115 VOLT, 60 HERTZ, 1 PHASE - 3.3 AMPS
- \*DOUBLE INSULATED POWER CORD WITH THREE PRONG GREEN DOT HOSPITAL PLUG.
- \*EXTERIOR 22 GAUGE STEEL POWDER COAT PAINT FINISH
- \*INTERIOR VACUUM FORM HIGH IMPACT POLYSTYRENE SEAMLESS WITH ROUNDED CORNERS FOR EASY CLEANING
- \*WARRANTY, 1 YEAR P&L, ADDITIONAL 4 YEARS ON COMPRESSORS PARTS ONLY

