

SO-LOW UNDERCOUNTER FREEZER MODEL MV25-4UCF

4.6 Cubic Ft. Upright – Temperature Range: 0°C to -25°C

One Solid Door

Specifications:

* Capacity	4.6 cubic ft.
* Temperature Range	0°C to -25°C
* Dimensions	
* Exterior w x d x h	24" x 24" x 33" 61 x 61 x 83 (cm)
* Interior w x d x h	20" x 19.5" x 24.25" 50 x 49 x 61 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Rear
Door	Single solid, locking
Shelves	Two
Insulation	Polyurethane foam
Refrigerant:	R404A
Compressor	1/3 hp
Defrost	Manual
Voltage: Domestic	115/60/1
Running Amps/Breaker	5/15
Supply Plug	5-15P NEMA
Shipping Weight	150 lbs.
Warranty-(Domestic)	One year parts and labor Compressors additional 4 years, no labor.

