

MODEL MV30-4UCF

LAB AND PHARMACY FREEZER -10°C TO -25°C

| | |
|---------------------------------|--------------------------------------|
| Capacity [Cu. Ft.] | 3.2 Cu. Ft. |
| Capacity [Liters] | 90 Liters |
| Temperature Range [°C] | -10°C to -25°C |
| Temperature Range [°F] | 14°F to -13°F |
| Dimensions | Width x Depth x Height |
| Exterior [Inches] | 27.75" x 27.75" x 34" |
| Interior [Inches] | 19.5" x 19.25" x 20.75" |
| Exterior [Centimeters] | 70 x 70 x 86 |
| Interior [Centimeters] | 50 x 49 x 53 |
| Compressor | Hermetic Compressor |
| Air Circulation | Direct Cooling |
| Defrost Method | Manual |
| Refrigerant | R600a |
| Exterior | Painted Steel |
| Interior | Painted Aluminum |
| Insulation | US EPA and SNAP approved |
| Lockable Door | Yes, Keyed |
| Porthole | Yes, 1" Diameter |
| Casters | 2, Rear |
| Leveling Legs | 2, Front |
| Storage | 2 Fixed Shelves, Base Shelf, Drawers |
| Shipping Weight | 165 lbs. |


Electrical

| | |
|-------------------------|---------------------------------|
| Voltage: | 115 Volts 60 Hertz 1 Phase |
| Amperage Line Running | 15 Amp Dedicated 1.66A |
| NEMA Configuration | NEMA 5-15 |
| Certification(s) | UL Listed Energy Star Listed |

Controller

| | |
|-----------------------|---|
| Display Type | Microprocessor, Digital |
| Alarms | High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure |
| Data Logging | Yes, Adjustable Intervals |
| Data Download | Yes, via USB, PDF Format |
| Min / Max Temperature | Yes, Display and Reset |
| Alarm Relay | Yes, Dry Contacts |
| Backup Battery | Yes, Only Powers Alarm |

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor