

SO-LOW UNDERCOUNTER REFRIGERATOR MODEL MV4-8UCRDA

8 Cubic Ft. (227 liters) Upright – Temperature Range: 1°C to 7°C / Solid Door

Specifications:

* Capacity	8 cubic ft. / 227 liters
* Temperature Range	1°C to 7°C
* Dimensions	
* Exterior w x d x h	29-7/8" x 24" x 34" 76 x 61 x 86 (cm)

Features:

Auto Defrost Cycle

Forced Air Ventilation through Front Grille

No Additional Ventilation Clearance Required
To Build in Under Counter

R-134A Refrigeration System

Keyed Door Lock

115 Volts, 60 Hertz, 1 Phase - 3.3 AMPS

Warranty, 1 Year Parts & Labor, Additional
4 Years on Compressors Parts Only

Double Insulated Power Cord With
Three Prong Green Dot Hospital
Plug

Exterior 22 Gauge Steel Powder Coat
Paint Finish

Interior Vacuum Form Impact
Polystyrene Seamless With Rounded
Corners For Easy Cleaning

CONTROL

Digital temperature control with two temperature displays. One display shows the set temperature, the other shows the actual temperature in the refrigerator.

ALARM

The battery operated alarm system will emit an audible and visual Signal when there is a mechanical or electrical failure. The alarm has An over and under temperature setting, alarm silencing switch, and Battery test switch.

REMOTE ALARM RELAY

A relay for a remote alarm hook-up is also provided.



Specifications subject to change without notice