

SO-LOW REFRIGERATED INCUBATOR

MODEL RIN10-20

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

20 Cubic Ft. Upright – Temperature Range: -10°C to 50°C / Solid Door

Specifications:

* Capacity	20 cubic ft.
* Temperature Range	-10°C to 50°C
* Dimensions	
* Exterior w x d x h	32" x 31" x 80.5" 81 x 78 x 204 (cm)
* Interior w x d x h	27" x 23.75" x 56.75" 68 x 60 x 144 (cm)

Features:

Temperature Control	Digital display with Hi/Low alarm
Doors	Single Solid, locking
Interior finish	White
Legs	Two leveling
Shelves	Five
Door Shelves	Six
Interior Lighting	Yes
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/4 hp
Thermostat	Adjustable
Defrost	Selectable on demand Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0/15
Supply Plug	5-15P NEMA
Shipping Weight	380 lbs.
Warranty-(Domestic)	1 ½ years parts and labor. Compressors additional 3 ½ years, no labor.

