

## *SO-LOW UPRIGHT FREEZER TO -80°C*

### *MODEL U80-30*

**30 Cubic Ft. (843 liters) Upright – Temperature Range: -40°C to -80°C**

#### **Upright Freezer**

##### ***Specifications:***

* Capacity	31 cubic ft. / 877 liters
* Temperature Range	-40°C to -80°C
* Dimensions	
* Exterior W x D x H	60.5" x 37" x 79.5" 153 x 93 x 201 (cm)
* Interior W x D x H	42" x 25" x 51.25" 106 x 63 x 130 (cm)

##### **TEMPERATURE CONTROL**

Microprocessor control displays set point and chamber temperature.

##### **ALARM SYSTEM WITH RELAY**

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

##### **REFRIGERATION SYSTEM**

A So-Low cascade system with two hermetic compressors. CFC & HCFC free refrigerants

##### **CONSTRUCTION**

Chamber is 14-ga. Zinc coated galvanized steel. Exterior is 16-ga. steel. Power coated cool gray.

##### **VOLTAGE**

20 AMP dedicated circuit required.  
ETL electrical approval.

##### **TECHNICAL SUPPORT**

Supplied by Nation-Wide-Network of service companies.

Specifications subject to change without notice

