

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SPECIAL PRICE: DEMO UNIT

SO-LOW LAB & PHARMACY GLASS DOOR FREEZER MODEL DH20-27GD-SA

27 Cubic Ft. Upright – Temperature Range: -15°C to -20°C / +5°F to -4°F

Specifications:

* Capacity	27 cubic ft.
* Temperature Range	-15°C to -20°C / +5°F to -4°F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	30" x 32" x 81" 76 x 81 x 206 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Doors	Single glass, locking
Interior finish	White powder coat, s/s floor
Casters	Four
Shelves	Five
Interior Lighting	One 20-watt fluorescent
Insulation	Urethane foam
Refrigerant:	R404A CFC free
Compressor	3/4 hp
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	13.8/20
Supply Plug	5-20P NEMA UL listed
Shipping Weight	450 lbs.
Warranty-(Domestic)	Three year parts and labor Compressor additional 2 years, no labor.



Specifications subject to change without notice