

## SPECIAL PRICE: DEMO UNIT

### SO-LOW LAB & PHARMACY FREEZER

### MODEL DH30-27SDSS-SA

**27 Cubic Ft. (764 liters) Upright – Temperature Range: -18°C to -25°C / 0°F to -13°F**

#### Specifications:

\* Capacity 27 cubic ft. / 764 liters  
\* Temperature Range -18°C to -25°C / 0°F to -13° F

**When ordering please specify Celsius or Fahrenheit Setting**

\* Dimensions  
\* Exterior w x d x h 30" x 33.75" x 81"  
76 x 85 x 206 (cm)

\* Interior w x d x h 27" x 28.5" x 55"  
69 x 72 x 140 (cm)

#### Features:

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Door	Single solid, locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four locking
Shelves	Three adjustable
Interior Lighting	One 20-watt fluorescent
Insulation	Urethane foam
Refrigerant:	R404A CFC free
Compressor	3/4 hp
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	13.8/20
Supply Plug	5-15P NEMA
Shipping Weight	440 lbs.
Warranty-(Domestic)	Three year parts and labor Compressor additional 2 years, no labor.

