

SPECIAL PRICE: OVERSTOCK

SO-LOW LAB & PHARMACY REFRIGERATOR

MODEL DH4-27GD-SA

27 Cubic Ft. (764 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

* Capacity	27 cubic ft. / 764 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	30" x 32" x 81" 76 x 81 x 206 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Door	Single glass
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four, locking
Shelves	Five, Adjustable
Interior Lighting	One
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	3.0 / 15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



NOTE:

UNIT DOES NOT
COME WITH
A LOCK

