

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SPECIAL PRICE: SCRATCH & DENT

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-38SGD-SA

38 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

* Capacity 38 cubic ft. / 1076 liters
* Temperature Range 2°C to 8°C / 36°F to 46°F

When ordering please specify Celsius or Fahrenheit Setting

* Dimensions
* Exterior w x d x h 43.5" x 32" x 81"
110 x 81 x 205 (cm)

* Interior w x d x h 40.5 x 27" x 55"
102 x 68 x 140 (cm)

Features:

Temperature Control	Digital display Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Doors	Double glass – sliding Lockable
Casters	Four, locking
Shelves	Ten, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane
Refrigerant:	R290A
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps	5.0/15
Supply Plug	5-15 NEMA UL listed
Shipping Weight	550 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice

