

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SPECIAL PRICE: SCRATCH & DENT

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-45SGD-SA

45 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

* Capacity 45 cubic ft. / 1286 liters

* Temperature Range 2°C to 8°C / 36°F to 46°F

When ordering please specify Celsius or Fahrenheit Setting

* Dimensions

* Exterior w x d x h 52" x 32" x 81"
132 x 81 x 206 (cm)

* Interior w x d x h 49" x 27" x 55"
125 x 69 x 140 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Doors	Two sliding glass, locking
Casters	Four, Locking
Shelves	Ten, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	610 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice