

## SPECIAL PRICE : OVERSTOCK

### SO-LOW LAB & PHARMACY FREEZER

### MODEL DHK30-23SD-SA

**23 Cubic Ft. Upright – Temperature Range: -15°C to -30°C**

#### Solid Door

#### Specifications:

|   |   |
|---|---|
| * Capacity  | 23 cubic ft.                            |
| * Temperature Range   | -15°C to -30°C                          |
| <b>*When ordering please specify Celsius or Fahrenheit Setting*</b> |   |
| * Dimensions  |   |
| * Exterior w x d x h  | 27" x 32" x 83.5"<br>68 x 81 x 212 (cm) |

#### Features:

|                      |                                      |
|----------------------|--------------------------------------|
| Temperature Control  | Digital display<br>With Hi/Low Alarm |
| Door                 | Single Solid, locking                |
| Access Port          | 1/2" Diameter                        |
| Interior Finish      | Stainless                            |
| Exterior Finish      | Stainless                            |
| Casters              | Four Casters                         |
| Shelves              | Three Adjustable                     |
| Interior Lighting    | LED Light                            |
| Insulation           | Polyurethane foam                    |
| Refrigerant:         | R600 CFC free                        |
| Air Circulation      | Forced air                           |
| Defrost              | Cycle, Auto                          |
| Voltage: Domestic    | 115/60/1                             |
| Running Amps/Breaker | 3 / 15                               |
| Supply Plug          | 5-15P NEMA<br>UL listed              |
| Shipping Weight      | 500 lbs.                             |



Specifications subject to change without notice