

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHK30-49SD-SA

49 Cubic Ft. Upright – Temperature Range: -15°C to -30°C

Solid Door

Specifications:

- * Capacity 49 cubic ft.
- * Temperature Range -15°C to -30°C
- *When ordering please specify Celsius or Fahrenheit Setting***

- * Dimensions
- * Exterior w x d x h 54" x 32" x 83.5"
137 x 81 x 212 (cm)

Features:

- | | |
|----------------------|--------------------------------------|
| Temperature Control | Digital display
With Hi/Low Alarm |
| Door | Single Solid, locking |
| Access Port | 1/2" Diameter |
| Interior Finish | Stainless |
| Exterior Finish | Stainless |
| Casters | Four Casters |
| Shelves | Six Adjustable |
| Interior Lighting | LED Light |
| Insulation | Polyurethane foam |
| Refrigerant: | R600 CFC free |
| Air Circulation | Forced air |
| Defrost | Cycle, Auto |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 9 / 15 |
| Supply Plug | 5-15P NEMA |
| | UL listed |
| Shipping Weight | 6500 lbs. |



Specifications subject to change without notice