

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## SPECIAL PRICE : OVERSTOCK

### *SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHK4-23GD-SA*

**23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F**

**Glass Door**

#### *Specifications:*

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
<b>*When ordering please specify Celsius or Fahrenheit Setting*</b>	
* Dimensions	
* Exterior w x d x h	27" x 32" x 83.5" 68 x 81 x 212 (cm)

#### *Features:*

Temperature Control	Digital display With Hi/Low Alarm
Door	Single Glass, locking
Access Port	1/2" Diameter
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four Casters
Shelves	Three Adjustable
Interior Lighting	LED Light
Insulation	Polyurethane foam
Refrigerant:	R600 CFC free
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	3 / 15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	500 lbs.



Specifications subject to change without notice