

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHK4-49GD-SA*

**49 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F**

**Glass Door**

### *Specifications:*

|   |  |
|---|--|
| * Capacity  | 49 cubic ft.                             |
| * Temperature Range   | 2°C to 8°C / 36°F to 46°F                |
| <b>*When ordering please specify Celsius or Fahrenheit Setting*</b> |  |
| * Dimensions  |  |
| * Exterior w x d x h  | 54" x 32" x 83.5"<br>137 x 81 x 212 (cm) |

### *Features:*

|                      |                                      |
|----------------------|--------------------------------------|
| Temperature Control  | Digital display<br>With Hi/Low Alarm |
| Doors                | Double Glass, locking                |
| Access Port          | 1/2" Diameter                        |
| Interior Finish      | Stainless                            |
| Exterior Finish      | Stainless                            |
| Casters              | Four Casters                         |
| Shelves              | Six Adjustable                       |
| Interior Lighting    | LED Light                            |
| Insulation           | Polyurethane foam                    |
| Refrigerant:         | R600 CFC free                        |
| Air Circulation      | Forced air                           |
| Defrost              | Cycle, Auto                          |
| Voltage: Domestic    | 115/60/1                             |
| Running Amps/Breaker | 6 / 15                               |
| Supply Plug          | 5-15P NEMA                           |
|                      | UL listed                            |
| Shipping Weight      | 650 lbs.                             |



Specifications subject to change without notice