

SPECIAL PRICE: SCRATCH & DENT

SO-LOW ECONOMY REFRIGERATOR

MODEL EC4-23GD-SA

23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C

Specifications:

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	27" x 32" x 78"
	68 x 81 x 198 (cm)
* Interior w x d x h	24" x 28.5" x 55"
	61 x 72 x 140 (cm)

Features:

Temperature Control	Digital display
Access Port	1/2" Diameter Top Back Wall
Doors	Single glass, locking
Interior finish	White powder coat
Legs	Four, leveling
Shelves	Five
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	3/15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	400 lbs.
Warranty-(Domestic)	Three year parts and labor
	Compressors additional
	2 years, no labor.

Specifications subject to change without notice

