

## SPECIAL PRICE - OVERSTOCK

### SO-LOW LAB & PHARMACY REFRIGERATOR

### MODEL DH4-12GD-SA

12 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

#### Glass Door

#### Specifications:

* Capacity	12 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
<b>*When ordering please specify Celsius or Fahrenheit Setting*</b>	
* Dimensions	
* Exterior w x d x h	25" x 23" x 66" 63 x 58 x 167 (cm)
* Interior w x d x h	22" x 18.5" x 44" 55 x 46 x 111 (cm)

#### Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Door	Single glass, locking
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Casters	Four locking
Shelves	Four Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	4.0 / 15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	400 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

