

SPECIAL PRICE – SCRATCH & DENT

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHK4-23SD-SA

23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	27" x 32" x 83.5" 68 x 81 x 212 (cm)
* Interior w x d x h	22.2" x 26.8" x 50.8" 56 x 68 x 129 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm
Door	Single Solid, locking
Access Port	1/2" Diameter
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four Casters
Shelves	Three Adjustable
Interior Lighting	LED Light
Insulation	Polyurethane foam
Refrigerant:	R600 CFC free
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	3 / 15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	500 lbs.
Warranty	2 Years Parts and Labor. Additional 3 years on compressor parts.



Specifications subject to change without notice

SPECIAL OFFER – SA03