

SPECIAL PRICE - OVERSTOCK

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHK4-23GD-SA

23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

- * Capacity 23 cubic ft.
- * Temperature Range 2°C to 8°C / 36°F to 46°F
- *When ordering please specify Celsius or Fahrenheit Setting***

* Dimensions

- * Exterior w x d x h 27" x 32" x 83.5"
68 x 81 x 212 (cm)
- * Interior w x d x h 22.2" x 26.8" x 50.8"
56 x 68 x 129 (cm)

Features:

- | | |
|----------------------|--|
| Temperature Control | Digital display
With Hi/Low Alarm |
| Door | Single Glass, locking |
| Access Port | 1/2" Diameter |
| Interior Finish | Stainless Steel |
| Exterior Finish | Stainless Steel |
| Casters | Four Casters |
| Shelves | Three Adjustable |
| Interior Lighting | LED Light |
| Insulation | Polyurethane foam |
| Refrigerant: | R600 CFC free |
| Air Circulation | Forced air |
| Defrost | Cycle, Auto |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 3 / 15 |
| Supply Plug | 5-15P NEMA
UL listed |
| Shipping Weight | 500 lbs. |
| Warranty | 2 Years Parts and Labor.
Additional 3 years on
compressor parts. |



Specifications subject to change without notice

SPECIAL OFFER # SA04