

SPECIAL PRICE – DEMO UNIT

SO-LOW LAB & PHARMACY FREEZER MODEL DH30-27SDSS

27 Cubic Ft. (764 liters) Upright – Temperature Range: -18°C to -25°C / 0°F to -13°F

Solid Door

Specifications:

* Capacity	27 cubic ft. / 764 liters
* Temperature Range	-18°C to -25°C / 0°F to -13°F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	30" x 33.75" x 81" 76 x 85 x 206 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

Features:

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Door	Single solid, locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four locking
Shelves	Three adjustable
Interior Lighting	One 20-watt fluorescent
Insulation	Urethane foam
Refrigerant:	R404A CFC free
Compressor	3/4 hp
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	13.8/20
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Three year parts and labor Compressor additional 2 years, no labor.



Specifications subject to change without notice

SPECIAL OFFER # SA06