

## SPECIAL PRICE – SCRATCH & DENT

### SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-49GDSS

49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Two Swinging Glass Doors

#### Specifications:

\* Capacity 49 cubic ft. / 1387 liters  
\* Temperature Range 2°C to 8°C / 36°F to 46°F  
**\*When ordering please specify Celsius or Fahrenheit Setting\***

\* Dimensions  
\* Exterior w x d x h 52" x 33.75" x 81"  
132 x 85 x 206 (cm)  
  
\* Interior w x d x h 49" x 28.5" x 55"  
125 x 72 x 140 (cm)

#### Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Doors	Two hinged glass, locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four, Locking
Shelves	Six, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	610 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice

## SPECIAL OFFER # SA11