

## SPECIAL PRICE – OVERSTOCK

### SO-LOW LAB & PHARMACY REFRIGERATOR

#### MODEL DHK4-49GD-SA

49 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Glass Door

#### Specifications:

- \* Capacity 49 cubic ft.
- \* Temperature Range 2°C to 8°C / 36°F to 46°F
- \*When ordering please specify Celsius or Fahrenheit Setting\*

#### \* Dimensions

- \* Exterior w x d x h 54" x 32" x 83.5"  
137 x 81 x 212 (cm)
- \* Interior w x d x h 49.4" x 26.8" x 50.8"  
125 x 68 x 129 (cm)

#### Features:

- Temperature Control Digital display  
With Hi/Low Alarm
- Doors Double Glass, locking
- Access Port 1/2" Diameter
- Interior Finish Stainless Steel
- Exterior Finish Stainless Steel
- Casters Four Casters
- Shelves Six Adjustable
- Interior Lighting LED Light
- Insulation Polyurethane foam
- Refrigerant: R600 CFC free
- Air Circulation Forced air
- Defrost Cycle, Auto
- Voltage: Domestic 115/60/1
- Running Amps/Breaker 6 / 15
- Supply Plug 5-15P NEMA  
UL listed
- Shipping Weight 650 lbs.
- Warranty 2 Years Parts and Labor.  
Additional 3 years on  
compressor parts.



Specifications subject to change without notice

## SPECIAL OFFER # SA15