

SPECIAL PRICE - OVERSTOCK

SO-LOW LAB & PHARMACY REFRIGERATOR

MODEL DHK4-49SD-SA

49 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Solid Door

Specifications:

* Capacity	49 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	54" x 32" x 83.5" 137 x 81 x 212 (cm)
* Interior w x d x h	49.4" x 26.8" x 50.8" 125 x 68 x 129 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm
Doors	Double Solid, locking
Access Port	1/2" Diameter
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four Casters
Shelves	Six Adjustable
Interior Lighting	LED Light
Insulation	Polyurethane foam
Refrigerant:	R600 CFC free
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	6 / 15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	650 lbs.
Warranty	2 Years Parts and Labor. Additional 3 years on compressor parts.



Specifications subject to change without notice

SPECIAL OFFER # SA16