

SPECIAL PRICE - SCRATCH & DENT

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-74GD-SA

72 Cubic Ft. (2095 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Specifications:

* Capacity	72 cubic ft. / 2095 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	75" x 32" x 81" 190 x 81 x 206 (cm)
* Interior w x d x h	72" x 27" x 61" 183 x 69 x 155 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Doors	Three glass – hinged – key lock
Shelves	Fifteen, adjustable
Casters	Four, Locking
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant	R290A CFC free
Compressor	1/3 hp
Air Circulation	Forced air
Defrost	Cycle Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	7.2 /15
Supply Plug	NEMA 5-15P UL listed
Shipping Weight	750 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice

SPECIAL OFFER # SA19