

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SPECIAL PRICE – SCRATCH & DENT

SO-LOW FLAMMABLE MATERIAL STORAGE FREEZER MODEL DHH20-30SDFMS - MANUAL DEFROST

Meets NFPA and OSHA guidelines for flammable material storage in general laboratory work areas

30 Cubic Ft. Upright to -20°C / Solid Door

Specifications:

* Capacity	30 cubic ft. / 849 liters
* Temperature	-15°C to -25°C
* Dimensions	
* Exterior w x d x h	35.75" x 35.5" x 72.5" 90.81 x 90.17 x 184.15 (cm)
* Interior w x d x h	32" x 22.5" x 65" 81 x 57 x 165 (cm)

Features:


*Thermostat	Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	4 fixed
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-404A
*Compressor	1/3 hp
*Air Circulation	Gravity
*Evaporator	Hot wall
*Defrost	Manual
*Voltage	115/60/1
*Running Amps	6
*Supply Plug	5-15 NEMA
*Shipping weight	378 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.



**IMAGES OF
DAMAGES
ON SECOND PAGE**

Specifications subject to change without notice

SPECIAL OFFER – SA # 22



ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS



SPECIAL OFFER – SA # 22