


MODEL CH40-15

15 CUBIC FOOT LABORATORY CHEST FREEZER


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	49.25 IN	21.63 IN	25.75 IN	125 CM	54 CM	65 CM
EXTERIOR DIMENSIONS	55.5 IN	29.75 IN	34.38 IN	140 CM	75 CM	87 CM
CRATED DIMENSIONS	64 IN	41 IN	56 IN	162 CM	104 CM	142 CM

NOTE THAT THE FREEZER LID CAN OPEN UP TO APPROXIMATELY 75° FROM HORIZONTAL POSITION

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 285 LBS.
SHIPPING WEIGHT: 325 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)
DOCK-TO-DOCK DELIVERY WITH LIFTGATE (ON PALLET)