


MODEL DH29-34SD

30 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT		WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	31 IN	26 IN	62 IN		78 CM	66 CM	157 CM
EXTERIOR DIMENSIONS	35 IN	33.75 IN	81 IN		89 CM	86 CM	209 CM
CRATED DIMENSIONS	46 IN	43 IN	86 IN		116 CM	109 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 400 LBS.

SHIPPING WEIGHT: 450 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)

DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)

INSIDE DELIVERY (WHITE GLOVE DELIVERY)