


MODEL DH4-38GD

38 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	43 IN	23.25 IN	60.5 IN	109 CM	59 CM	153 CM
EXTERIOR DIMENSIONS	47.25 IN	29.25 IN	78.5 IN	120 CM	74 CM	199 CM
CRATED DIMENSIONS	53 IN	43 IN	86 IN	134 CM	109 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 450 LBS.
SHIPPING WEIGHT: 500 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)
DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)
INSIDE DELIVERY (WHITE GLOVE DELIVERY)