


MODEL DH4-49GD-SS

49 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


| | WIDTH | DEPTH | HEIGHT | | WIDTH | DEPTH | HEIGHT |
|---------------------|---------|----------|--------|--|--------|--------|--------|
| INTERIOR DIMENSIONS | 47.5 IN | 28.5 IN | 55 IN | | 120 CM | 72 CM | 140 CM |
| EXTERIOR DIMENSIONS | 52 IN | 33.75 IN | 81 IN | | 132 CM | 85 CM | 206 CM |
| CRATED DIMENSIONS | 56 IN | 43 IN | 86 IN | | 142 CM | 109 CM | 218 CM |

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 575 LBS.

SHIPPING WEIGHT: 610 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)

DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)

INSIDE DELIVERY (WHITE GLOVE DELIVERY)