


MODEL DH4-49SD-SS

49 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	47.5 IN	28.5 IN	55 IN	120 CM	72 CM	140 CM
EXTERIOR DIMENSIONS	52 IN	33.75 IN	81 IN	132 CM	85 CM	206 CM
CRATED DIMENSIONS	56 IN	43 IN	86 IN	142 CM	109 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 575 LBS.
SHIPPING WEIGHT: 610 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)
 DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)
 INSIDE DELIVERY (WHITE GLOVE DELIVERY)