


# MODEL DHK4-10GD

## 11.12 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


|                     | WIDTH    | DEPTH   | HEIGHT  | WIDTH | DEPTH | HEIGHT |
|---------------------|----------|---------|---------|-------|-------|--------|
| INTERIOR DIMENSIONS | 22.75 IN | 21 IN   | 44 IN   | 58 CM | 53 CM | 112 CM |
| EXTERIOR DIMENSIONS | 25.75 IN | 26.5 IN | 69.5 IN | 65 CM | 67 CM | 176 CM |
| CRATED DIMENSIONS   | 32 IN    | 32 IN   | 80 IN   | 81 CM | 81 CM | 203 CM |

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

### SIDE VIEW



### FRONT VIEW



# SO-LOW

AUTHORIZED DEALER OF **LIEBHERR**

[WWW.SO-LOW.COM](http://WWW.SO-LOW.COM)



**NET WEIGHT : 190 LBS.**

**SHIPPING WEIGHT: 240 LBS.**

### SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )

DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )

INSIDE DELIVERY ( WHITE GLOVE DELIVERY )