


MODEL DHK4-49GD-T

46.62 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT		WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	52.75 IN	26.25 IN	59.5 IN		134 CM	67 CM	152 CM
EXTERIOR DIMENSIONS	57.25 IN	34.75 IN	78.75 IN		145 CM	88 CM	200 CM
CRATED DIMENSIONS	85 IN	41 IN	86 IN		215 CM	104 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 750 LBS.
SHIPPING WEIGHT: 800 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)
DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)
INSIDE DELIVERY (WHITE GLOVE DELIVERY)