


MODEL LRT50G2HC

46.4 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	48.5 IN	25.5 IN	59 IN	48.5 CM	25.5 CM	59 CM
EXTERIOR DIMENSIONS	56.5 IN	32.75 IN	85.25 IN	143 CM	83 CM	216 CM
CRATED DIMENSIONS	85 IN	41 IN	86 IN	215 CM	104 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

AUTHORIZED DELAER OF **LIEBHERR**

WWW.SO-LOW.COM



NET WEIGHT : 485 LBS.

SHIPPING WEIGHT: 580 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)

DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)

INSIDE DELIVERY (WHITE GLOVE DELIVERY)