


# MODEL U85-13

## 13 CUBIC FOOT ULTRA LOW TEMPERATURE FREEZER


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	20 IN	22 IN	51.25 IN	50 CM	55 CM	130 CM
EXTERIOR DIMENSIONS	35 IN	34.5 IN	79.5 IN	88 CM	87 CM	201 CM
CRATED DIMENSIONS	46 IN	43 IN	86 IN	116 CM	109 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

**SIDE VIEW**



**FRONT VIEW**



\*TO PASS THROUGH DOORWAYS, THE DOOR CAN OPEN A FULL 180°, TO REDUCE THE UNIT EXTERIOR DEPTH TO 30.5"

**SO-LOW**

[WWW.SO-LOW.COM](http://WWW.SO-LOW.COM)



**NET WEIGHT : 675 LBS.**

**SHIPPING WEIGHT: 850 LBS.**

**SHIPPING OPTIONS**

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )

DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )

INSIDE DELIVERY ( WHITE GLOVE DELIVERY )