

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

CHEST FREEZER RACKS

FOR ALL STANDARD 3.00-INCH STORAGE BOXES

The highest quality racks, constructed from corrosion resistant stainless steel, including handles.
Rack listed is for standard 3.00-inch high cardboard, aluminum stainless steel, or plastic boxes.
Custom racks available upon request.

CHEST STYLE RACK – 13-4-3

SPECIFICATIONS:

Rack Dimensions	Inches	Centimeters
Width	5.5"	14
Depth	5.5"	14
Height	12.75"	32

Box Dimensions

Width	5.25"	13
Depth	5.25"	13
Height	3"	7

Box Capacity: 4 Boxes (Included)
Box Cell Divider 81-Cell Divider (Included)
Rack Construction: Stainless Steel
Box Construction: Cardboard

Maximum Rack Capacity

2 Cu. Ft. Chest Freezer can hold up to 6 Racks



If you have any questions, or need any additional information, please contact us via phone at 513-772-9410 or via email at sales@so-low.com. Please visit our website at www.so-low.com for additional details regarding these products.