

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## CHEST FREEZER RACKS

### FOR ALL STANDARD 2.00-INCH STORAGE BOXES

The highest quality racks, constructed from corrosion resistant stainless steel, including handles.  
Rack listed is for standard 2.00-inch high cardboard, aluminum stainless steel, or plastic boxes.  
Custom racks available upon request.

## CHEST STYLE RACK – 13-6-2

### SPECIFICATIONS:

<b>Rack Dimensions</b>	Inches	Centimeters
Width	5.75"	14
Depth	5.5"	14
Height	13.25"	33

### Box Dimensions

Width	5.25"	13
Depth	5.25"	13
Height	2"	5

Box Capacity: 6 Boxes ( Included )  
Box Cell Divider 81-Cell Divider ( Included )  
Rack Construction: Stainless Steel  
Box Construction: Cardboard

### Maximum Rack Capacity

2 Cu. Ft. Chest Freezer can hold up to 6 Racks



If you have any questions, or need any additional information, please contact us via phone at 513-772-9410 or via email at [sales@so-low.com](mailto:sales@so-low.com). Please visit our website at [www.so-low.com](http://www.so-low.com) for additional details regarding these products.