

CHEST FREEZER RACKS

FOR ALL STANDARD 3.00-INCH STORAGE BOXES

For All Models of Chest Freezers and Refrigerators

The highest quality racks, constructed from corrosion resistant stainless steel, including handles. Racks listed are for standard 3.00-inch high cardboard, aluminum stainless steel, or plastic boxes. Custom racks available upon request.

CHEST STYLE RACK – 16-5-3

Specifications:

Rack Dimensions

- Width (A): 5.75"
- Height (B): 16.0"
- Depth (C): 5.5"
- Box Capacity: 5 Boxes

Rack includes boxes and 81 cell divider.

Box Dimensions

- Width: 5.25"
- Depth: 5.25"
- Height: 3.0"

Rack Capacity:

- 3 cu. ft. chest can hold 8 racks
- 5 cu. ft. chest can hold 15 racks
- 9 cu. ft. chest can hold 16 racks
- 12 cu. ft. chest can hold 27 racks
- 14 cu. ft. chest can hold 30 racks
- 17 cu. ft. chest can hold 36 racks
- 22 cu. ft. chest can hold 36 racks
- 27 cu. ft. chest can hold 48 racks

