

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

UPRIGHT FREEZER RACKS

FOR ALL STANDARD 2.00-INCH STORAGE BOXES

The highest quality racks, constructed from corrosion resistant stainless steel, including handles.

Rack listed is for standard 2.00-inch high cardboard, aluminum stainless steel, or plastic boxes.

Custom racks available upon request.

UPRIGHT DRAWER STYLE RACK – 17-12-2-DR

SPECIFICATIONS:

Rack Dimensions	Inches	Centimeters
Width	16.5"	41
Depth	5.5"	14
Height	9.5"	24

Box Dimensions

Width	5.25"	13
Depth	5.25"	13
Height	2"	5

Box Capacity: 12 Boxes (Included)
Box Cell Divider 81-Cell Divider (Included)
Rack Construction: Stainless Steel
Box Construction: Cardboard

Maximum Rack Capacity

13 Cu. Ft. Upright Freezer can hold up to 15 Racks
18 Cu. Ft. Upright Freezer can hold up to 20 Racks
22 Cu. Ft. Upright Freezer can hold up to 25 Racks
25 Cu. Ft. Upright Freezer can hold up to 30 Racks
28 Cu. Ft. Upright Freezer can hold up to 30 Racks
30 Cu. Ft. Upright Freezer can hold up to 35 Racks



If you have any questions, or need any additional information, please contact us via phone at 513-772-9410 or via email at sales@so-low.com. Please visit our website at www.so-low.com for additional details regarding these products.