

CHEST FREEZER RACKS

FOR ALL STANDARD 2.00-INCH STORAGE BOXES

The highest quality racks, constructed from corrosion resistant stainless steel, including handles.
Rack listed is for standard 2.00-inch high cardboard, aluminum stainless steel, or plastic boxes.
Custom racks available upon request.

CHEST STYLE RACK – 29-13-2

SPECIFICATIONS:

Rack Dimensions	Inches	Centimeters
Width	5.75"	14
Depth	5.5"	14
Height	28.25"	71


Box Dimensions

Width	5.25"	13
Depth	5.25"	13
Height	2"	5

Box Capacity: 13 Boxes (Included)
Box Cell Divider 81-Cell Divider (Included)
Rack Construction: Stainless Steel
Box Construction: Cardboard

Maximum Rack Capacity

3 Cu. Ft. Chest Freezer can hold up to 8 Racks
5 Cu. Ft. Chest Freezer can hold up to 15 Racks
9 Cu. Ft. Chest Freezer can hold up to 16 Racks
12 Cu. Ft. Chest Freezer can hold up to 27 Racks
14 Cu. Ft. Chest Freezer can hold up to 30 Racks
17 Cu. Ft. Chest Freezer can hold up to 36 Racks
22 Cu. Ft. Chest Freezer can hold up to 36 Racks
27 Cu. Ft. Chest Freezer can hold up to 48 Racks



If you have any questions, or need any additional information, please contact us via phone at 513-772-9410 or via email at sales@so-low.com. Please visit our website at www.so-low.com for additional details regarding these products.