

MODEL C40-21

CHEST ULTRA LOW FREEZER 0°C TO -40°C

Capacity [Cu. Ft.] 21 Cu. Ft.
Capacity [Liters] 594 Liters
Temperature Range [°C] 0°C to -40°C
Temperature Range [°F] 32°F to -40°F

Dimensions Width x Depth x Height
 Exterior [Inches] 91" x 36" x 49"
 Interior [Inches] 81" x 21.5" x 21"
 Exterior [Centimeters] 231 x 91 x 124
 Interior [Centimeters] 205 x 54 x 53

Refrigeration System Single Stage System
 Compressors 1, Hermetic Compressor
 Air Circulation Gravity Flow
 Defrost Method Manual
 Refrigerant HCFC & CFC Free
 Exterior Powder Coated Steel
 Interior Zinc Coated Galvanized Steel
 Insulation Polyurethane (Ecomate) Foam
 Additional Insulation Beaded Foam Sub-Lids
 Lockable Lid Lockable with Padlock
 Porthole Optional 1/2"
 Casters 8, Swivel with Locks
 Storage Base Shelf
 Shipping Weight 1025 Lbs.

Electrical
 Certification(s) ETL Listed

Controller
 Display Type Microprocessor, Digital
 Alarms High/Low Temp, Power Failure, Low Battery
 Alarm Relay Yes, Dry Contacts
 Automatic Restart Yes, Automatically Restarts After Power Failure
 Backup Battery Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor