

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW CHEST STYLE FREEZER TO -40°C

MODEL C40-5

5 Cubic Ft. (142 liters) Chest – Temperature Range: 0°C to -40°C
Top Opening

Specifications:

* Capacity	5 cubic ft. / 142 liters
* Temperature Range	0°C to -40°C
* Dimensions	
* Exterior W x D x H	40.5" x 31" x 47" 102 x 81 x 119 (cm)
* Interior W x D x H	30" x 18" x 16" 76 x 46 x 40 (cm)

TEMPERATURE CONTROL

Microprocessor control displays set point and chamber temperature.

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low single stage system with one hermetic compressor. CFC & HCFC free refrigerants.

CONSTRUCTION

Chamber is 14-ga. Zinc coated galvanized steel. Exterior is 16-ga. steel. Power coated cool gray.

VOLTAGE

15 AMP dedicated circuit required.
ETL electrical approval.

TECHNICAL SUPPORT

Supplied by Nation-Wide-Network of service companies.



Specifications subject to change without notice