

SO-LOW CHEST STYLE FREEZER TO -80°C

MODEL C80-27

27 Cubic Ft. (766 liters) Chest – Temperature Range: -40°C to -85°C

Top Opening

Specifications:

* Capacity	27 cubic ft. / 766 liters
* Temperature Range	-40°C to -80°C
* Dimensions	
* Exterior W x D x H	98" x 36" x 43.5" 248 x 91 x 110 (cm)
* Interior W x D x H	67" x 23" x 31.25" 170 x 57 x 79 (cm)

TEMPERATURE CONTROL

Microprocessor control displays set point and chamber temperature.

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low cascade system with two hermetic compressors. CFC & HCFC free refrigerants

CONSTRUCTION

Chamber is 14-ga. Zinc coated galvanized steel. Exterior is 16-ga. steel. Power coated cool gray.

VOLTAGE

208 volts 60 hertz, 1 phase.

20 AMP dedicated circuit required.

ETL electrical approval.

TECHNICAL SUPPORT

Supplied by Nation-Wide-Network of service companies.

Specifications subject to change without notice

