

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW CHEST STYLE FREEZER TO -25°C MODEL CH25-15*

**15 Cubic Ft. (424 liters) Chest – Temperature Range: 0°C to -25°C  
Top Opening**

### ***Specifications:***

*Capacity	15 cubic ft. / 424 liters
*Temperature Range	0°C to -25°C
*Dimensions	
- Exterior w x d x h	55.50" x 29.75" x 34.38"
- Interior w x d x h	49.25" x 21.63" x 25.75"

### **CONTROL**

Digital Temperature control with two Temperature displays. One display shows the set temperature, the other shows the actual temperature in the freezer.

### **ALARM**

Battery operated alarm system provides An audible and visual signal in a high/low temperature alarm condition. The alarm will also signal in the event of a power failure. Also provided is an alarm silencing and battery test feature.

MANUAL DEFROST

CASTERS FOR EASY MOBILITY

1/3HP HERMETIC MOTOR COMPRESSOR

115 VOLT, 60 HERTZ, 1 PHASE

WARRANTY One Year Parts  
90 Days Labor

### **TECHNICAL SUPPORT**

Supplied by Nation-Wide-Network of service companies.

**ETL Electrical Approval**

