

SO-LOW MODEL CH25-20

CHEST STYLE FREEZER TO -25°C

Capacity [Cu. Ft.]	20 Cubic Feet
Capacity [Liters]	566 Liters
Temperature Range [°C]	-0°C to -25°C
Temperature Range [°F]	+32°F to -13°F

Dimensions (W x D x H)

Exterior [Inches] 74" x 30" x 34"

Exterior [Centimeters] 188 x 76 x 86

Interior [Inches] 67" x 21.5" x 25"

Interior [Centimeters] 170 x 54 x 63

Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage: 15 Amp Dedicated Line

2 Running Amps

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

Construction:

Temperature Control Digital Display
High / Low Alarm
Min / Max Temperature

Alarm Relay Dry Contacts

Defrost Manual

Opening Top Opening Lid

Lockable Lid/Door No

Interior / Exterior Steel

Insulation Polyurethane Foam

Compressors Hermetic Compressor

Primary Refrigerant(s) R134A

Air Circulation Gravity Flow

Casters 4, Lockable

Shelving None

Approximate BTU/Hour 600

Approximate Decibels 55

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & 90 Days Labor